





ABOUT US



CONTENTS

4-5



COLLAR

7-10



SHOULDER & ELBOW

12-17



WRIST

19-25



CORSET

27-30



KNEE & HIP

32-36




ANKLE

38-39



FOOT HEALTH

41-42



COMPRESSION STOCKING

44



OTHER PRODUCTS

46



PILLOW

48



WALKING AID



NAVIGATING YOUR PATH TO OPTIMAL HEALTH

COLLAR







101 SOFT FOAM COLLAR

It is produced by cotton sheathed sponge fabric in accordance with the anatomical neck structure. With its elastic structure, it restricts neck movements, controls them and keeps the soft tissues loose by keeping the area warm.

Indications

- Used in soft tissue crushes
- Used in neck eclipses
- Used At the beginning of neck hernia
- Used for thermotherapy



S	M	L	XL
28-33	34-38	39-43	44-50
PERIMETER/ CM			
COLOR  BLACK  BEIGE			



102 Cervical Collar With Chins

It is made of plastic material, designed according to the anatomical shape of the neck. Limits neck movements with its chin support. Height adjustment can be made with the help of the Velcro straps on the top.

Indications

- Used in neck hernia
- It is used for post-operative protection and stabilization



S	M	L	XL	XXL
32-35	36-39	40-43	43-46	46-50
PERIMETER/ CM				



103 Cervical Collar Without Chins

It is made of plastic material, designed according to the anatomical shape of the neck. Height adjustment can be made with the help of the Velcro straps on the top. It reduces pressure on the neck muscles by 5\MRK the RecO area.

Indications

- Used in neck hernia
- It is used for post-operative protection and stabilization



S	M	L	XL	XXL
32-35	36-39	40-43	43-46	46-50
PERIMETER/ CM				



104 NELSON COLLAR

It is made of durable plastozote material with a thickness of 12 mm in accordance with the anatomical neck structure. It controls neck movements, especially by restricting it to the front and back, and keeps the soft tissues loose by keeping the area warm.

Indications

- Used in soft tissue crushes
- Used in neck eclipses
- Used in neck hernia
- Used for post-traumatic protection
- Used for thermotherapy



S	M	L	XL
28-33	34-38	39-43	44-50
PERIMETER/ CM			

COLLAR



105 PHILADELPHIA COLLAR

It is produced from 12 mm thick plastozote material in accordance with the anatomical neck structure. The neck collar, which is applied to the neck in the form of two parts, completely restricts neck movements and keeps the soft tissues loose by keeping the area warm.

Indications

- Used in soft tissue crushes
- Used in vertebral fractures
- Used in neck hernia
- First aid to trauma
- Used in post-surgical rehabilitation



S	M	L
33-37	38-42	42-46
PERIMETER/ CM		



106 Philadelphia Tracheostomy Collar

It is produced from 12 mm thick plastozote material in accordance with the anatomical neck structure. The neck collar, which is applied to the neck in the form of two parts, completely restricts neck movements and keeps the soft tissues loose by keeping the area warm.

Indications

- Used in soft tissue crushes
- Used in vertebral fractures
- Used in neck hernia
- First aid to trauma
- Used in post-surgical rehabilitation



S	M	L
33-37	38-42	42-46
PERIMETER/ CM		



108 MINERVA ORTHOSIS

It is made of polyethylene plastic and cotton fabric material. The neck collar, which is applied to the neck and upper body in a multi-part manner, restricts the neck movements in the desired position by taking support from the shoulder and upper body elements. With the help of its adjustable mechanism, it helps

Indications

- Used in neck hernia
- Used for post-operative protection and stabilization



S	M	L
28-33	34-39	40-45
PERIMETER/ CM		



111 Foam Collar with Plastic Support

It is produced by cotton sheathed sponge fabric in accordance with the anatomical neck structure. With its

Indications

- Used in soft tissue crushes
- Used in neck eclipses
- Used at the beginning of neck hernia
- Used for thermotherapy



S	M	L	XL
28-33	34-38	39-43	44-50
PERIMETER/ CM			

SHOULDER & ELBOW



SHOULDER & ELBOW



152 arm sling (mesh)

It enables the patient to keep his/her arm close to the body and parallel to the ground. Suitable for left and right arm use.

Indications

- Used in shoulder dislocations and semi-dislocations
- Treatment of forearm diseases
- Used for resting after surgery
- It is used to prevent edema formation in the fracture



CODE	S	M	L	XL	XXL
152	30-33	33-36	36-39	39-42	42-45
153	ONE SIZE				
PERIMETER/CM					



155 arm sling

Its inner surface is made of terry cloth and its outer surface is sponge fabric and it is designed in accordance with arm and shoulder anatomy. It ensures that the patient's arm is stable, parallel to the ground, and close to the body. Suitable for left and right arm use.

Indications

- Used in shoulder dislocations and semi-dislocations
- Treatment of forearm diseases
- Used for resting after surgery
- It is used to prevent edema formation in the fracture



S	M	L	XL	XXL
30-33	33-36	36-39	39-42	42-45
PERIMETER/CM				



156 arm sling (simple)

Its inner surface is made of terry cloth and its outer surface is sponge fabric and it is designed in accordance with arm and shoulder anatomy. With the help of its practical design, it allows the patient to keep his arm close to the body and parallel to the ground. Suitable for left and right arm use.

Indications

- Used in shoulder dislocations and semi-dislocations
- Treatment of forearm diseases
- Used for resting after surgery
- It is used to prevent edema formation in the fracture

**ONE
SIZE**



157 arm sling (30°)

It is a padded arm sling at a 30 degree angle. It is made of sponge fabric and with the help of its design, it allows the patient's arm to stay parallel to the ground close to the body. Suitable for left and right arm use. There is an additional belt for the waist circumference to facilitate daily use.

Indications

- Used in rotator cuff muscle tears
- Used for rehabilitation after shoulder implant prosthesis
- Used for rehabilitation after frozen shoulder arthroscopy

**ONE
SIZE**



158 Arm Sling (45°-60°)

It is a padded arm sling at a 45-60 degree angle. It is made of sponge fabric and with the help of its design, it allows the patient's arm to stay parallel to the ground close to the body. Suitable for left and right arm use. There is an additional belt for the waist circumference to facilitate daily use.

Indications

- Used in rotator cuff muscle tears
- Used for rehabilitation after shoulder implant prosthesis
- Used for rehabilitation after frozen shoulder arthroscopy

**ONE
SIZE**



159 Clavicle Bandage

It is made of sponge fabric with a Terry cloth on the inner surface. With the Velcro adjustment, it allows the shoulder to be pulled back with the desired tension and provides a smooth posture. With the help of its clavicle bandage design, it provides the necessary immobilization and prevents dislocation.

Indications

- Used in Clavicle (collarbone) fractures
- Used in dislocations of the joint between the clavicle (collarbone) and sternum (breastbone)
- Used in shoulder joint injuries and postural disorders (postural kyphosis)



CODE	XS	S	M	L	XL	XXL
159	18-24	25-30	31-36	37-42	43-48	49-54
160	ONE SIZE					
PERIMETER/CM						



161 VELPEAU BANDAGE

It is made of sponge fabric and is designed in accordance with the arm and shoulder anatomy. It provides immobilization of the shoulder and upper arm region in abduction. Bandage wrapping around the chest area enhances the effect of the Velcro on the joint. The ease of use is a plus.

Indications

- Used in Shoulder dislocations and semi-dislocations
- Used in tears of the rotator cuff muscles
- Used in Fractures of the humerus (arm bone) and scapula
- Used for the soft tissue injuries in the shoulder and arm



CODE	S	M	L	XL	XXL
161	30-33	33-36	36-39	39-42	42-45
162	ONE SIZE				
PERIMETER/CM					



165 shoulder support

It is made of neoprene material and designed in accordance with the anatomy of the shoulder. It is a shoulder bandage that protects the heat in the area and increases blood circulation thanks to the compression it applies to the shoulder area. Thanks to the applied pressure, it prevents edema. Suitable for use in the following cases:

Indications

- Used in Shoulder bone fractures
- Used in Shoulder arthritics
- Used in post-operative irritations
- Used for rehabilitation after surgery

**ONE
SIZE**

SHOULDER & ELBOW




167 HUMERUS BRACE

It is a shoulder orthosis made of sponge fabric with a Terry cloth on the inner surface. The plastic and aluminum parts protect the humerus in shoulder dislocations, after surgery and plastercast. An optional arm hanger is available.

Indications

- Used for protection and stabilization in shoulder dislocations
- Used in humeral head fractures
- Used for rehabilitation after surgery
- Used for rehabilitation after shoulder plaster application



		L R	
S	M	L	
22-24	24-27	27-30	
PERIMETER/ CM			




168 Contractor Arm Orthosis

It is made of thermoplastic material and is designed in accordance with shoulder and arm anatomy. With the help of the mechanism in the elbow area, the desired adjustment can be made. It can be adjusted in degrees of extension and flexion. It prevents the wrist from falling, length adjustment can be made on the orthosis according to the arm length. Shoulder strap provides comfort in use.

Indications

- Used After ligament injuries
- Used After elbow implant prosthesis application
- Used to speed up the healing process
- Used for conservative treatment after elbow dislocation



		L R	
S	M	L	
22-24	24-27	27-30	
PERIMETER/ CM			



169 epicondylitis elbow support

It is made of neoprene material and is designed in accordance with the anatomy of the elbow. With the help of the special 3-point silicone pad that provides local compression, it reduces the pain in the tendon. Designed for left and right arm use, it accelerates blood circulation and preserves body temperature.

Indications

- Used for Lateral epicondylitis
- Used in painful conditions of the elbow



CODE	S	M	L	XL	XXL
169	22-24	24-26	26-28	28-30	30-32
170	ONE SIZE				
PERIMETER/CM					



174 elbow immobilizer

It is made of breathable neoprene fabric. It prevents sweating thanks to towel laminated inside surface. Thanks to the steel bars, it ensures that the arm stops in stabilization and the immobilization of the elbow.

Indications

- Elbow joint traumas.
- Rest after surgery.
- It helps protection and immobilization.

**ONE
SIZE**





175 Epicondylitis Bandage

It is made of neoprene fabric and designed in accordance with the anatomy of the elbow. Thanks to its special pad that provides local compression, it reduces the pain in the tendon. Designed for left and right arm use, it accelerates blood circulation and preserves body temperature. It is Universal size and can be adjusted individually.

Indications

- Tendinitis
- Golfer's elbow
- Tennis elbow
- It is used for post-operative protection and rest.

**ONE
SIZE**



226 KNITTED ELBOW SUPPORT

- It is made of neoprene fabric and designed in accordance with the anatomy of the elbow. It reduces tendon tension and pain with local compression.

Indications

- Provides protection and support to the elbow
- Arthritis
- It is suitable for daily use against pain



S	M	L	XL	XXL
22-24	24-26	26-28	28-30	30-32
PERIMETER/ CM				



227 Knitted Epicondylitis Elbow Support

- It is made of neoprene fabric and designed in accordance with the anatomy of the elbow. It provides compression on the muscle groups attached to the lateral epicondyle.

Indications

- In lateral epicondylitis
- It is used in painful conditions of the elbow



S	M	L	XL	XXL
22-24	24-26	26-28	28-30	30-32
PERIMETER/ CM				





indications

-



Indications

-



Indications

- | | | |
|---|-------|-------|
|  | S-M | L-XL |
| | 15-19 | 19-25 |
| PERIMETER / CM | | |



Indications

- | | S-M | L-XL |
|----------------|-------|-------|
| | 15-19 | 19-25 |
| PERIMETER / CM | | |

HAND WRIST



204 Wrist Splint

It is made of breathable anti-allergic fabric and is designed in accordance with the anatomy of the wrist. It has an anatomically shaped aluminum support that provides immobilization of the wrist in neutral position. There is a bandage that increases the stabilization around the wrist.

Indications

- Carpal tunnel syndrome
- Tendinitis, tenosynovitis, tenovaginitis
- Used for resting in painful conditions of the wrist
- Used for injury and post-operative protection



**ONE
SIZE**



205 Wrist Splint with Thumb

It is made of breathable anti-allergic fabric and is designed in accordance with the anatomy of the wrist. It has an anatomically shaped aluminum support that provides immobilization of the wrist in neutral position. There is a bandage that increases the stabilization around the wrist.

Indications

- Carpal tunnel syndrome
- Tendinitis, tenosynovitis, tenovaginitis
- Used for resting in painful conditions of the wrist
- Used for injury and post-operative protection



**ONE
SIZE**



206 Lace Wrist Splint

-t Ms QaHe oJ aMr te\JaFrMc [Mth FreathaFle Qesh. -t TroZMHes MQQoFMIM^atMoR oJ the [rMst MR ReYtral TosMtMoR. 8here Ms aR e\tra FaRHake oR the [rMst that MRcreases staFMIM^atMoR. -t Ms laceH Jor ease oJ Yse.

Indications

- Carpal tunnel syndrome
- Tendinitis, tenosynovitis, tenovaginitis
- Used for resting in painful conditions of the wrist
- Used for injury and post-operative protection



**ONE
SIZE**



207 Lace Wrist Splint With Thumb

-t Ms QaHe oJ aMr te\JaFrMc [Mth FreathaFle Qesh. -t TroZMHes MQQoFMIM^atMoR oJ the [rMst aRH thYQF MR ReYtral TosMtMoR. 8here Ms aR e\tra FaRHake oR the [rMst that MRcreases staFMIM^atMoR. -t Ms laceH Jor ease oJ Yse.

Indications

- Used When carpal tunnel and dequervain are together
- g 9seH to +raHe NoMRt calcM5catMoRs
- g 9seH Jor 4rotectMoR MR thYQF IMKaQeRt la\Mt]
- Used after plaster in fractures



**ONE
SIZE**





208 Knitted Wrist Splint

-t Ms TroHYceH YsMRK se\MFle FreathaFle (ORMttMRK techRoloK]. 8he TlastMc sYTTort aroYRH the YIRar steloidreducespressureonvesselsandnerves while protectingthe bonyprominence. There issupport onthevolarfacetthat provides immobilization of the wrist in neutral position. The bandage is designed for Qa\MQYQ 5t coQJor t aRH eJlectMZe Yse QaMRtaMRs FoH] teQTeratYre accelerates FlooH cMrcYlatMoR aRH healing process.

Indications

- Carpal tunnel syndrome
- Tendinitis, tenosynovitis, tenovaginitis
- Painful conditions of the wrist
- Injuryand post-surgery



S	M	L	XL	XXL
15-17	17-19	19-21	21-23	23-25
PERIMETER/ CM				



210 DeQuervain Splint

Itismadeofbreathable anti-allergic fabric. It is designed in accordance with the anatomy of the hand, wristandthumb. Withits elasticstructure, it strengthens the musclesby compressingthe tissuesand TreZeRts s[ellMRK. ;Mth the heLT oJ aLYQMRYQ sTIMRt Mt 5\es the thYQF aRH restrMcts Mts QoZeQeRt.

Indications

- Used in Dequervain syndrome

ONE
SIZE



211 Wrist Bandage

Itismadeofbreathable anti-allergic fabric. It is designed in accordance with the anatomy of the hand, wristandthumb. Withits elasticstructure, it strengthens the musclesby compressingthe tissuesand TreZeRts s[ellMRK. ;Mth the heLT oJ aLYQMRYQ sTIMRt Mt 5\es the thYQF aRH restrMcts Mts QoZeQeRt.

Indications

- Used in Dequervain syndrome

ONE
SIZE



211-1 Wrist Bandage Woven Fabric

It is made of breathable anti-allergic woven fabric. It is designed in accordance with the anatomy of the hand,wrist andthumb. Withitselasticstructure,it strengthens the musclesby compressingthe tissues aRH TreZeRts s[ellMRK. ;Mth the heLT oJ aLYQMRYQ sTIMRt Mt 5\es the thYQF aRH restrMcts Mts QoZeQeRt.

Indications

- Used in Dequervain syndrome

ONE
SIZE

COLOR  BLACK  BEIGE

HAND WRIST



215 The Frog Splint

7TecMal aLYQMRYQ 5RKer sTIMRts [Mth Tlasto^ote coateH MRRer sYrJace. -t TroZMHes MQQoFMIM^atMoR oj the 5RKer MR e\teRsMoR. -t Ms YseH Jor 5RKer JractYres teRHoR aRH RerZe MRNYrMes straMRs aRH sTraMRs.

Indications

- Used for Finger fractures
- Used for Tendon and nerve injuries
- Used for Finger injuries and sprains



S	M	L
4,5-5	5-6	6-7
PERIMETER/ CM		



216 The Baseball Splint

7TecMal aLYQMRYQ 5RKer sTIMRts [Mth Tlasto^ote coateH MRRer sYrJace. -t TroZMHes MQQoFMIM^atMoR oj the 5RKer MR e\teRsMoR. -t Ms YseH Jor 5RKer JractYres teRHoR aRH RerZe MRNYrMes MRNYrMes aRH sTraMRs.

Indications

- Used for Finger fractures
- Used for Tendon and nerve injuries
- Used for Finger injuries and sprains



S	M	L
4,5-5	5-6	6-7
PERIMETER/ CM		



217 buttonier splint

-t Ms TroHYceH JroQ therQoTlastMc QaterMal aRH HesMKReH MR accorHaRce [Mth the 5RKer aRatoQ]. -t Ms a 5RKer or thosMs that OeeTs the 5RKer MR e\teRsMoR accorHMRK to the ToMRt TrMRcMTle.

Indications

- Used for Rehabilitation after surgery



S	M	L	XL
4	5	5-6	6-7
PERIMETER / CM			



218 Mallet Finger Splint

-t Ms QaHe oj sTecMal TlastMc QaterMal aRH Ms HesMKReH MR accorHaRce [Mth 5RKer aRatoQ]. -t Ms a 5RKer or thosMs that OeeTs the 5RKer NoMRt MR e\teRsMoR. 8he 5RKer Ms Jree JroQ the 4-4 NoMRt. -t Ms YseH MR e\teRsor tendon injuries, fractures, sprains and sprains.

Indications

- Used in PIP joint fractures, strains, and sprains



1	2	3	4	5	6
4,5-5	5-5,5	5,5-6	6-6,5	6,5-7	7-8
PERIMETER/CM					




219 Thermoplastic Splint


It is made of thermoplastic material and its inner surface is covered with terry cloth. Keeps the wrist and 5Rkers MR the restMRK TosMtMoR. -t Ms YseH to staFMIM^e the haRH MR the restMRK TosMtMoR after a Tlaster cast or in painful conditions of the wrist.

Indications

- Used After plaster cast in bone fractures joining the wrist
- Used For resting soft tissue injuries
- Carpal Tunnel Syndrome
- Used for Rehabilitation after surgery
- Used after Tendon and soft tissue damage



S	M	L
14-18	18-22	22-26
PERIMETER/ CM		





221 Thermoplastic Antispasticity Splint

It is made of thermoplastic material and its inner surface is covered with terry cloth. It is a hand and wrist orthosis made of thermoplastic material, the inner surface of which is covered with plastozote. Due to its corrYKateH 5Rker HesMKR aRH resMstaRce to teRsMoR Mt Ms esTecMall] YseH MR cases oJ sTastMcMt] MR the haRH.

Indications

- Spasticity
- Wrist discomforts



S	M	L
14-18	18-22	22-26
PERIMETER/ CM		





222 Thermoplastic Dynamic Splint


It is made of thermoplastic material and its inner surface is covered with terry cloth. With the help of its haRH te\tYre Mt MQQoFMIM^es the [rMst MR the aRatoQMcal TosMtMoR. ;Mth the heIT oJ aLYQMRYQ aRH sTrMRK aTTaratYs oTeR 5Rkers Ho FeRHMRK e\ercMses.

Indications

- Used for Correcting muscle inactivity
- Used for correcting joint instability
- Used In ulnar nerve damage
- Used in post-surgical rehabilitation



S	M	L
14-18	18-22	22-26
PERIMETER/ CM		





223 Thermoplastic Cleinert Splint


It is made of thermoplastic material and its inner surface is covered with terry cloth. The splint, designed MR accorHaRce [Mth the aRatoQ] oJ the [rMst allo[s TassMZe se\MoR aRH actMZe e\teRsMoR. -t allo[s earl] coRtrollH QoZeQeRt MR se\or teRHoR MRNYrMes aRH ajter teRHoR reTaMr.

Indications

- Used for Rehabilitation after surgery



S	M	L
14-18	18-22	22-26
PERIMETER/ CM		



HAND WRIST



224 Middle Finger Fracture Splint

4roZMHes staFMIM^atMoR MR the [rMst aRH QMHHe 5RKer. -t caR Fe aHNYsteH to the HesMreH tMKhtRes [Mth the Zelcro straTs oR Mt. 8here Ms aR aLYQMRYQ sTIMRt MRsMHe Jor Qa\MQYQ staFMIM^atMoR aRH MQQoFMIM^atMoR.

Indications

- Used for Finger fractures
- Used for Tendon and nerve injuries
- Used for Finger injuries and sprains
- Used for Rehabilitation after surgery

**ONE
SIZE**



225 Boxer Fracture Splint

-t TroZMHes staFMIM^atMoR MR the [rMst IMttle 5RKer aRH rMRK 5RKers. -t caR Fe aHNYsteH to the HesMreH tMKhtRes [Mth the Zelcro straTs oR Mt. 8here are t[o aLYQMRYQ sTIMRts MRsMHe Jor Qa\MQYQ staFMIM^atMoR aRH MQQoFM-lization.

Indications

- MP joint fractures
- Used for Rehabilitation after surgery

**ONE
SIZE**

CORSET



CORSET



250 CHEST CORSET

-t Ms QaHe oj se\MFle FreathaFle aRtM TersTMrMRK aRtM allerKM Mc 5shMRK IMRe JaFrMc. 8he chest corset sYTTorts thesternumareawith terrycloth. Thereare twounderwiresonthechest.Itsupportsandprotectsthe sternumbyreducing pain caused by sudden movement.

Indications

- Used In rib fractures
- Used for upper body stabilization
- Used after thoracic surgeries and traumas



CODE	S	M	L	XL	XXL
250	75-85	85-95	95-105	105-115	115-125
251	ONE SIZE				
PERIMETER/CM					



252 Umbilical Hernia Corset

-t Ms QaHe oj se\MFle FreathaFle aRtM TersTMrMRK cottoR JaFrMc. -t has a sTecMal HesMKReH TaH that TroZMHes compression to the herniated area. It is applied in umbilical hernia.

Indications

- Used In abdominal soft tissue hernias
- Used for organ failure
- Supportingthe abdominal wall
- Used in post-surgical rehabilitation



CODE	S	M	L	XL	XXL
252	75-85	85-95	95-105	105-115	115-125
253	ONE SIZE				
PERIMETER/CM					



254 Colostomy Corset (26cm)

-t Ms QaHe oj se\MFle FreathaFle aRtM TersTMrMRK aRtM allerKM c)lastMc IMRe JaFrMc. -t Ms YseH to reHYce the riskofherniaaftersurgeryincolostomycases.Itsupportsthecolostomybagandalsotheabdomen.Itis produced with 26 cm.

Indications

- Protecting the colostomy bag
- Providing support after ostomy surgery
- Used toreduce the risk of parastomal hernia after surgery



S	M	L	XL	XXL
75-85	85-95	95-105	105-115	115-125
PERIMETER / CM				



255 Colostomy Corset (16cm)

-t Ms QaHe oj se\MFle FreathaFle aRtM TersTMrMRK aRtM allerKM c)lastMc IMRe JaFrMc. -t Ms YseH to reHYce the riskofherniaaftersurgeryincolostomycases.Itsupportsthecolostomybagandalsotheabdomen.Itis produced with 16 cm.

Indications

- Protecting the colostomy bag
- Providing support after ostomy surgery
- Used to reduce the risk of parastomal hernia after surgery



S	M	L	XL	XXL
75-85	85-95	95-105	105-115	115-125
PERIMETER / CM				



256 Abdominal Corset (16cm)

-t Ms QaHe oj seVMFle FreathaFle aRtM TersTMrMRK aRtM allerKM Mc elastMc IMRe JaFrMc ;Mth Mts elastMc strYctYre Mt protects and strengthens the muscles by compressing the tissues in the abdomen, and prevents swelling. It is produced with 16 cm.

Indications

- Used As a postpartum recovery
- Used for organ failure
- Supporting the abdominal wall
- Used in post-surgical rehabilitation



CODE	S	M	L	XL	XXL
256	75-85	85-95	95-105	105-115	115-125
257	ONE SIZE				
PERIMETER/CM					



258 Abdominal Corset (26cm)


-t Ms QaHe oj seVMFle FreathaFle aRtM TersTMrMRK aRtM allerKM Mc elastMc IMRe JaFrMc ;Mth Mts elastMc strYctYre Mt protects and strengthens the muscles by compressing the tissues in the abdomen, and prevents swelling. It is produced with 26 cm.

Indications

- Used As a postpartum recovery
- Used for organ failure
- Supporting the abdominal wall
- Used in post-surgical rehabilitation



CODE	S	M	L	XL	XXL
258	75-85	85-95	95-105	105-115	115-125
259	ONE SIZE				
PERIMETER/CM					



260 ELASTIC CORSET

With its elastic structure, it compresses the tissues, strengthens the muscles and maintains body temperature. It can be used underneath daily clothing.

Indications

- Used In abdominal organ cold
- Used for low back pain
- Used for thermotherapy



S	M	L	XL	XXL
75-85	85-95	95-105	105-115	115-125
PERIMETER / CM				



261 ELASTIC CORSET WITH HOOK

With its elastic structure, it compresses the tissues, strengthens the muscles and maintains body temperature. It can be used underneath daily clothing.

Indications

- Used In abdominal organ cold
- Used for low back pain
- Used for thermotherapy



S	M	L	XL	XXL
75-85	85-95	95-105	105-115	115-125
PERIMETER / CM				

CORSET



262 POSTUREX CORSET

-t Ms QaHe o] se\MFle FreathaFle RoR s[eatMRK aRtM allerKMc elastMc IMRe JaFrMc. -t streRktheRs the QYscles byprovidingcompression to thetissuesinthe backandwaist area with itselastic structure. Withthehelp ofplasticunderwires on the back, it restricts the movement of the back vertebrae and puts it at rests.

Indications

- On the hunchback
- Used for Back muscle weakness
- Used in posture disorder



CODE	S	M	L	XL	XXL
262	75-85	85-95	95-105	105-115	115-125
263	ONE SIZE				
PERIMETER/CM					



264 MATERNITY CORSET

Itismadeofbreathable and anti-perspiring elastic cotton material. Used before and after birth. This postpartumcorsethelps to recoverlargesagging in the anterior abdominal region. It gives a feeling of coR5HeRce aRH coQJor t Jor QoZeQeRt MR HaMI] IMJe.

Indications

- Helps Back and waist pain during pregnancy
- Baby at risk of miscarriage
- Used in organ failure



CODE	S	M	L
264	70-80	80-90	90-100
265	ONE SIZE		
PERIMETER/ CM			



266 Lumbar Corset with Belt (26cm)

-t Ms QaHe o] se\MFle FreathaFle aRtM allerKMc elastMc cQ elastMc IMRe JaFrMc. ;Mth Mts elastMc strYctYre Mt strengthensthemusclesby compressing thetissues in the waistandhipsandpreventsswelling. Withthe helpof4steelunder wires, it restricts the movement of the lumbar vertebrae and supports the vertebrae andputsthebackat rest.

Indications

- Used In lumbar spine hernias
- Used in Sudden and long-term low back pain
- Used in post-surgical rehabilitation



CODE	S	M	L	XL	XXL
266	75-85	85-95	95-105	105-115	115-125
267	ONE SIZE				
PERIMETER/CM					



269 Lumbar Corset with Belt (32cm)

-t Ms QaHe o] se\MFle FreathaFle aRtM allerKMc elastMc cQ elastMc IMRe JaFrMc. ;Mth Mts elastMc strYctYre Mt strengthensthemusclesby compressing thetissues in the waistandhipsandpreventsswelling. Withthe helpof4steelunder wires, it restricts the movement of the lumbar vertebrae and supports the vertebrae andputsthebackat rest.

Indications

- Used In lumbar spine hernias
- Used in Sudden and long-term low back pain
- Used in post-surgical rehabilitation



CODE	S	M	L	XL	XXL
269	75-85	85-95	95-105	105-115	115-125
270	ONE SIZE				
PERIMETER/CM					



271 DORSOLUMBAR CORSET

-t Ms TroHYceH JroQ se\MFle FreathaFle aRtM allerKMc elastMc IMRe JaFrMc. ;Mth Mts elastMc strYctYre Mt streRKth-ensthemusclesbycompressing thetissues inthe back,waistand hips.With the steelunder wires on bothsidesofthespine, it helps to carry the upper body and rests by restricting the movement of the lumbar vertebrae appropriately.

Indications

- Used for spinal traumas
- Used in osteoporosis
- Used in Sudden and long-term back and waist pain
- Used in post-surgical rehabilitation



CODE	S	M	L	XL	XXL
271	75-85	85-95	95-105	105-115	115-125
272	ONE SIZE				
PERIMETER/CM					



273 Neoprene Abdominal Corset (26cm)

-t Ms QaHe oJ se\MFle FreathaFle ReoTreRe JaFrMc. ;Mth Mts elastMc strYctYre Mt Troprotects aRH streRKtheRs the muscles by compressing the tissues in the abdomen, and prevents swelling. It is produced with 26 cm.

Indications

- Used As a postpartum recovery
- Used for organ failure
- Supporting the abdominal wall
- Used in post-surgical rehabilitation

**ONE
SIZE**



274 T.L.S.O. CORSET

Itismadeofsponge fabric, the inner part of which is terry cloth, and there is a thermoplastic material thatisportableand provides rigid support to the abdomen and hard part of the corset. The support in the sternum,whichisdifferent in this corset, which can beattached and removed likethe Dorso-lumbar corset TroZMHes se\MoR e\teRsMoR lateral teRHeRc] JYll coRtrol aRH MQQoFMIM^atMoR oJ the trYRO MR the Tost oT TerMoH. -t caR Fe e\teRHeH or shor teReH accorHMRK to the FoH] leRKth at the sterRYQ.

Indications

- Used After severe trauma
- Used In vertebral fractures
- Onthe hunchback
- Used in post-surgical rehabilitation



S	M	L
60-85	85-110	110-140
PERIMETER/ CM		



275 Hyperextension Corset (JEweet)

It is a body orthosis that is made of special aluminum, the surfaces in contact with the skin are covered [Mth terr] cloth aRH accorHMRK to the ToMRt TrMRcMTle Mt restrMcts the se\MoR oJ the FoH] aRH eRcoYraKes Mt MR the h]Tere\teRsMoR TosMitMoR. -t IMQMts the QoZeQeRts oJ the sTMRe to a certaMR e\teRt [Mth the heIT oJ Mts screwsystemonthe orthosis,lengthandwidthadjustmentscanbemade.

Indications

- Used for Fractures and traumas in the spine
- Used for kyphosis
- Used in post-operative rehabilitation



S	M	L
75-95	95-115	115-125
PERIMETER/ CM		

CORSET



278 Hernia Truss (Double)


It is designed according to the anatomical shape of the thigh region. It is made of tough elastic cotton fabric. There is a special support that applies pressure to the groin area. It is used in cases of unilateral inguinal hernia.

Indications

- Used in inguinal hernia



S	M	L
70-90	90-110	110-130
PERIMETER/ CM		



279 Hernia Truss (SINGLE)


It is designed according to the anatomical shape of the thigh region. It is made of tough elastic cotton fabric. There is a special support that applies pressure to the groin area. It is used in cases of unilateral inguinal hernia.

Indications

- Used in inguinal hernia



LEFT & RIGHT		
S	M	L
70-90	90-110	110-130
PERIMETER/ CM		



280 POSTUREX CORSET COMFORT

It provides compression to the tissues in the back and waist area with its elastic structure. With the help of plastic underwires on the back, it restricts the movement of the back vertebrae and puts it at rest.

Indications

- Used for Back muscle weakness
- Used in posture disorder

**ONE
SIZE**



282 Lumbosacral Corset With Belt Comfort (26cm)

It is manufactured from breathable fabric. The height of corset is 26 cm. It supports waist and until to low-back. There are two plastic supports at the back of corset and two portable steel supports both sides of spine. 2 crossbelts from the back to the navel area, provide additional stabilisation to the lumbar area.

Indications

- In sacralization and lumbalization
- Sacroiliac joint problems
- Sudden and prolonged low back pain
- Rehabilitation after surgical intervention



S-M	L-XL
75-100	100-125
PERIMETER/ CM	





283 Abdominal Corset 3 Bands

-t Ms QaHe oj se\MFle FreathaFle aRtM TersTMrMRK aRtM allerKM elastMc IMRe JaFrMc. ;Mth Mts elastMc strYctYre Mt protectsandstrengthensthemuscles bycompressing the tissuesintheabdomen,andpreventsswelling. It has 3 bands, all bands can apply the force in the region of use independently.

Indications

- Used As a postpartum recovery
- Used for organ failure
- Supporting the abdominal wall
- Used in post-surgical rehabilitation



S-M	L-XL
75-100	100-120
PERIMETER / CM	
COLOR  BLACK  BEIGE	





284 Abdominal Corset 4 Bands

-t Ms QaHe oj se\MFle FreathaFle aRtM TersTMrMRK aRtM allerKM elastMc IMRe JaFrMc. ;Mth Mts elastMc strYctYre Mt protectsandstrengthensthemuscles bycompressing the tissuesintheabdomen,andpreventsswelling. It has 4 bands, all bands can apply the force in the region of use independently.

Indications

- Used As a postpartum recovery
- Used for organ failure
- Supporting the abdominal wall
- Used in post-surgical rehabilitation



S-M	L-XL
75-100	100-120
PERIMETER / CM	
COLOR  BLACK  BEIGE	



285 Knitted Lumbosacral Corset

-t Ms TroHYceH YsMRK se\MFle FreathaFle (ORMttMRK techRoloK) aRH Ms HesMKReH MR accorHaRce [Mth the waistanatomy. There are 4 steel underwires at the back. Steel underwires can be tilted manually according to the lumbar lordosis. It adapts perfectly to the body with its conformity to the anatomical structure and superior quality. It reduces pain by reducing the tension in the muscles in the lumbar region.

Indications

- In sacralization and lumbarization
- In sacroiliac joint problems
- Sudden and long-term low back pain
- In rehabilitation after surgery



S	M	L	XL	XXL
75-85	85-95	95-105	105-115	115-125
PERIMETER / CM				



286 Knitted Lumbosacral Corset with Pad

-t Ms TroHYceH YsMRK se\MFle FreathaFle (ORMttMRK techRoloK) aRH Ms HesMKReH MR accorHaRce [Mth the waistanatomy. It is a corset with a specially designed triangular pad on the back lumbar part that creates a massage effect. Increases circulation with massage effect. It stabilizes the sacroiliac joint and the lumbosacral joint, and the compression applied reduces

Indications

- In sacralization and lumbarization
- In sacroiliac joint problems
- Sudden and long-term low back pain
- In rehabilitation after surgery



S	M	L	XL	XXL
75-85	85-95	95-105	105-115	115-125
PERIMETER / CM				

CORSET



290 ELASTIC TRICOT CORSET

With its elastic structure, it compresses the tissues, strengthens the muscles and maintains body temperature. It can be used underneath daily clothing.

Indications

- Used in abdominal organ cold
- Used for low back pain
- Used for thermotherapy



S	M	L	XL	XXL
75-85	85-95	95-105	105-115	115-125
PERIMETER / CM				



291 Hip Corset

It is made of neoprene fabric. It is used as a support for warmth and compression in muscle tensions around the thigh and hip.

Indications

- Muscle strains around the thigh and hip
- Used as support for warmth and compression



S	M	L	XL	XXL
75-85	85-95	95-105	105-115	115-125
PERIMETER / CM				
COLOR <input checked="" type="radio"/> BLACK <input type="radio"/> BEIGE				



281 WORKING CORSET (16CM)

It is 16cm high corset manufactured from air fabric and supporting to waist and back. Consist 4 plastic supports at the back. 2 cross belts from the back to the navel area provides stabilization.

Indications

- In sacralization and lumbalization
- Sacroiliac joint problems
- Sudden and prolonged low back pain
- Rehabilitation after surgical intervention



S-M	L-XL-XXL
70-90	90-120
PERIMETER / CM	

KNEE & HIP



KNEE & HIP



300 Elastic Knee Support

It is made of acrylic fabric. With its elastic structure, it compresses the tissues, strengthens the muscles and maintains body temperature.

Indications

- Used for knee pain
- Used for thermotherapy



S	M	L	XL	XXL
31-34	34-37	37-40	40-44	44-48
PERIMETER/CM				



301 Knee Support Closed

It is made of neoprene fabric. It strengthens the muscles by compressing the tissues with its elastic structure and prevents swelling. It restricts the closing and opening movements of the knee and protects it from possible damage.

Indications

- Used for Post-operative protection and rehabilitation
- Used For protection during sports activities
- Provides heat and compression



CODE	S	M	L	XL	XXL
301	31-34	34-37	37-40	40-44	44-48
302	ONE SIZE				
PERIMETER/CM					



305 Knee Support with Patella

It is made of neoprene fabric. It is designed to support the knee joint. Patella is open around and supported by a silicone pad. With the help of the heat it spreads in the area and the compression it provides, it reduces the pain. It is suitable for use before and after surgery. It can be preferred for daily use as it does not put pressure on the kneecap.

Indications

- Used for Kneecap pain
- Used for Post-operative protection and rehabilitation
- Used For protection during sports activities
- Provides heat and compression



CODE	S	M	L	XL	XXL
305	31-34	34-37	37-40	40-44	44-48
306	ONE SIZE				
PERIMETER/CM					



307 Knee Support with Patella & Ligament

It is made of neoprene fabric. It is designed to support the knee joint. Patella circumference is open and TaHHeH. *le\MFle YRHer [Mres oR the sMHe MRcrease staFMIM^atMoR aRH relMeZe the loaH oR the ORee. -t JYrther increases the stabilization with the Velcro tape on the thigh.

Indications

- Arthritis, bursitis, tendinitis
- Used In sprains and strains
- Used for meniscus problems
- Used for Ligament Instability
- For post-operative protection



CODE	S	M	L	XL	XXL
307	31-34	34-37	37-40	40-44	44-48
308	ONE SIZE				
PERIMETER/CM					



309 Knee Support with Patella and Hinges

It is made of neoprene material. Patella circumference is open and supported by a silicone pad. There are steel NoMRts or Foth sMHes TroZMHMRK Qa\MQYQ sYTTort aRH allo[MRK coRtrolleH QoZeQeRt oJ the ORee. -t TreZeRts ORee h]Tere\teRsMoR. -t has :elcro straTs that MRcrease staFMIM^atMoR oR the thMKh aRH calj.

Indications

- Used for ligament injuries and ligament instability
- Used for meniscus problems
- Used In sprains and strains
- Restricts knee movement after surgery and speeds up the healing process



CODE	S	M	L	XL	XXL
309	31-34	34-37	37-40	40-44	44-48
310	ONE SIZE				
PERIMETER/CM					



311 Knee Support with Crossed Ligaments

It is made of neoprene fabric. Designed to support knee joints, Patella circumference is open and support- eH [Mth a sMIMcoRe TaH. ;Mth the heIT oJ seQM se\MFle sYTTor ts oR Mts MRRer aRH oYter sMHes Mt reImeZes the loaH oR the lateral IMKaQeRts. -t has aR e\tra sYTTortMRK FaRHKe JroQ the oYtsMHe to the MRsMHe [Mth the 3-point principle. It strengthens the muscles by compressing the tissues with its elastic structure, provides stabilization in the anterior cruciate ligaments and patella.

Indications

- Used in Anterior cruciate ligament injuries and knee instability
- Used for Patellar stabilization
- Used for meniscus problems
- Used In ligament damage
- Used for Post-operative protection



CODE	S	M	L	XL	XXL
311	31-34	34-37	37-40	40-44	44-48
326	ONE SIZE				
PERIMETER/CM					



314 Patellar Band

It is made of neoprene material. It has an anatomically designed pad that supports the patellar tendon. It does not interfere with joint movements. The compression effect it applies reduces the tension and pain in the tendon. It can be used comfortably all day. It can be applied to both strings.

Indications

- Patellar tendinitis
- Os-good Schallter disease
- Used for Patella femoral pain
- Chondromalacia patella

**ONE
SIZE**



315 Knitted Knee Brace

-t Ms TroHYceH YsMRK se\MFle FreathaFle (ORMttMRK techRoloK] aRH Ms HesMKReH MR accorHaRce [Mth the knee anatomy. The heat effect it provides on the knee joint provides long-term relief on the aching joints, muscles and surrounding tissues. Daily as it does not put pressure on the kneecap preferable for use.

Indications

- Knee cap pain
- Meniscus problems
- First degree sprains and ligament instability
- Post-operative protection



S	M	L	XL	XXL
31-34	34-37	37-40	40-44	44-48
PERIMETER/CM				

KNEE & HIP



316 Knitted Knee Support With Patella

-t Ms TroHYceH YsMRK 5e\MFle FreathaFle (ORMttMRK techRoloK) aRH Ms HesMKReH MR accorHaRce [Mth the kneeanatomy. The silicone pad around the patella relieves the pressure on the patella. Thanks to the heat it spreads in the area and the compression it provides, it reduces the pain. Ideal for people with kneecap pain when climbing, climbing or squatting stairs.

Indications

- RcreaseH MRtra artMcYlar 5YMH
- In rehabilitation after surgery
- Protection during and after sports activities
- Kneepain
- Arthritis



S	M	L	XL	XXL
31-34	34-37	37-40	40-44	44-48
PERIMETER/ CM				



317 Knitted Knee Support With Patella & Ligament

-t Ms TroHYceH YsMRK 5e\MFle FreathaFle (ORMttMRK techRoloK) aRH Ms HesMKReH MR accorHaRce [Mth the ORee aRatoQ]. 4atella Ms sYTTorteH F] a sMIMcoRe TaH. *le\MFle YRHer [Mres oR Foth sMHes MRcrease staFMIM - zation and thus relieve the load on the knee. Thanks to the heat and compression it spreads in the area, it reduces the pain. It is suitable for use before and after surgery.

Indications

- RcreaseH MRtra artMcYlar 5YMH
- In rehabilitation after surgery
- In lateral ligament damage and kneecap pain
- Meniscus problems
- Arthritis, bursitis



S	M	L	XL	XXL
31-34	34-37	37-40	40-44	44-48
PERIMETER/ CM				



318 Knee Immobilizer

It is made of sponge fabric with a terry cloth on the inner surface. It provides immobilization of the knee in the R s MoR. -t has reQoZaFle steel sYTTor ts aRH :elcro that TroZMHe TosterMor QeHMal aRH lateral staFMIM^atMoR. Supportive for knee injuries and post-surgery, with the help of its three-piece pad design, it can be easily adjusted for leg circumference. It is available in 45 cm, 55 cm, and 65 cm lengths and is standard size.

Indications

- Used for Post-surgery in meniscus injuries
- Used for lateral ligament damage
- Used for Post-operative protection and immobilization
- Used for Rehabilitation after plaster cast



S	M	L
45	55	65
PERIMETER/ CM		



319 Knee Brace with Adjustable Angle

It is an adjustable knee orthosis with an aluminum bar designed in accordance with the anatomical leg shape aRH terr] cloth sYrJaces MR coRtact [Mth the sOMR. *le\MoR q q q q q q q q q q aRH e\teRsMoR q q q q q q q q q . -t TroZMHes e\celleRt staFMIM^atMoR aRH sYTTort after knee injuries.

Indications

- Controlled post-operative injuries in ligament injuries and graded joint mobilization
- Post-surgery for meniscus injuries
- At rest at the desired angle after the surgical operation
- Used In knee recurvatum injuries



ADJUSTABLE	
SHORT	LONG
47	58
PERIMETER/ CM	



321 Thigh Support

It is made of neoprene fabric. Designed to support the thigh, with the help of the heat and compression it spreads to the area, it reduces the pain.

Indications

- Used in muscle strains
- Used in painful conditions of the thigh
- For protection during and after sports activities
- Used to provide warmth and compression

ONE
SIZE



323 Calf Support

It is made of neoprene fabric. Designed to support the calf, with the help of the heat and compression it spreads to the area, it reduces the pain.

Indications

- Used in muscle strains
- Used in painful conditions of the calf
- For protection during and after sports activities
- Used to provide warmth and compression

ONE
SIZE



370 Hip Brace

It is a hip abduction orthosis that provides stabilization of the hip joint made of anatomically designed thermoplastic material. The inner surface in contact with the skin is covered with terry cloth. It can be used after total hip replacement or revision, for stabilization purposes, and for recurrent hip dislocation cases.

Indications

- Used after total hip replacement or revision
- Used for stabilization purposes
- Used for recurrent hip dislocation cases

ONE
SIZE



327 KNEE SUPPORT ELASTIC LONG

It is made of terry fabric. With its elastic structure, it compresses the tissues, strengthens the muscles and maintains body temperature.

Indications

- Used for knee pain
- Used for thermotherapy



S	M	L	XL	XXL
31-34	34-37	37-40	40-44	44-48
PERIMETER/ CM				

ANKLE





350 Eight Bandage

It is made of neoprene fabric and is designed in accordance with the foot anatomy. It relieves the foot structure, relieves the edema and reduces the pain.

Indications

- Ankle sprains
- To support ligament weakness
- Mild soft tissue bruises
- Used for protection in sports activities

ONE
SIZE



351 ANKLE SUPPORT

It is made of neoprene fabric and is designed in accordance with the foot anatomy. It relieves the foot structure, relieves the edema and reduces the pain.

Indications

- Mild ligament injuries and instability
- For protective and support purposes after surgery and injury
- It is used in cases of edema, hematoma and effusion



CODE	S	M	L	XL	XXL
351	21-23	23-25	25-27	27-29	29-31
352	ONE SIZE				
PERIMETER/CM					



353 Ankle Support with Malleolar Pad

It is made of neoprene fabric and is designed in accordance with the foot anatomy. There are pads on both sides of the ankle that cross the inner and outer malleolar. It relieves the foot structure, relieves the edema and reduces the pain.

Indications

- Mild ligament injuries and instability
- For protective and support purposes after surgery and injury
- It is used in cases of edema, hematoma and effusion



CODE	S	M	L	XL	XXL
353	21-23	23-25	25-27	27-29	29-31
354	ONE SIZE				
PERIMETER/CM					



355 Ankle Support with Ligament

It is made of neoprene fabric and is designed in accordance with the foot anatomy. Anatomically designed with an elastic crossband, which is wrapped in the shape of eight. With the compression effect it applies, it supports the ankle and reduces edema and pain.

Indications

- In ligament instability and tears
- In ankle sprains
- In post-plaster recovery
- Used for protection or support after surgery and injury



CODE	S	M	L	XL	XXL
355	21-23	23-25	25-27	27-29	29-31
356	ONE SIZE				
PERIMETER/CM					

ANKLE



357 Ankle Support with Plastic Stay

Designed for ankle injury, rehabilitation and sports medicine. Plastic stay on both sides increases ankle stability. Reinforced crisscross strap pulls stay toward ankle to prevent sprain or swelling. Open heel applies less pressure on the Achilles tendon area.

Indications

- In ligament instability and tears
- In ankle sprains
- In post-plaster recovery
- Used for protection or support after surgery and injury



CODE	S	M	L	XL	XXL
357	21-23	23-25	25-27	27-29	29-31
358	ONE SIZE				
PERIMETER/CM					



359 Knitted Ankle Support

*le\MFle FreathaFle TroHYceH YsMRK (ORMttMRK techRoloK) aRH HesMKReH MR accorHaRce [Mth the Joot anatomy. The heat effect it provides on the ankle provides long-term relief on aching joints, muscles and surrounding tissues. The compression effect it applies reduces edema and pain.

Indications

- In ankle sprains and strains
- To provide warmth and compression
- For compression purposes in post-injury edema



S	M	L	XL	XXL
21-23	23-25	25-27	27-29	29-31
PERIMETER/CM				



360 Knitted Ankle Support With Malleolar Pad

*le\MFle FreathaFle TroHYceH YsMRK (ORMttMRK techRoloK) aRH HesMKReH MR accorHaRce [Mth the Joot anatomy. There are silicone pads around the malleoli. By supporting the ankle, it reduces pain and edema with its compression effect. It is comfortable to use.

Indications

- In ankle sprains
- In rehabilitation after surgery
- In rehabilitation after plaster cast
- In ligament damage
- For compression purposes in post-injury edema



S	M	L	XL	XXL
21-23	23-25	25-27	27-29	29-31
PERIMETER/CM				



361 Knitted Ankle Support With Achilles Tendon

*le\MFle FreathaFle TroHYceH YsMRK (ORMttMRK techRoloK) aRH HesMKReH MR accorHaRce [Mth the Joot anatomy. There is silicone support on both sides of the Achilles tendon and under the heel. It also has a portable heel placed under the other heel to prevent imbalance between both legs. It protects the heel and Achilles tendon, reduces the tension in the tendon.

Indications

- Achilles tendon pain
 - After surgery or injury
 - Retrocalcaneal bursitis
- g %cYte aRH chroRMc MR5aQQatMoRs of the %chMlles teRHOR



S	M	L	XL	XXL
21-23	23-25	25-27	27-29	29-31
PERIMETER/CM				



362 DORSI-FLEXION BANDAGE

It is made of neoprene fabric and is designed in accordance with the foot anatomy. It is an ankle brace that improves the patient's gait.

Indications

- Low foot deformity
- Peroneal neuropathy
- Muscle weakness in paralyzed patients
- Rehabilitation after surgery

ADJUSTABLE



363 Aircast Ankle Support

The outer part is made of special thermoplastic material, and the inner part has pads with air cushions.

Indications

- In recurrent ankle sprains in ligament damage
- Rehabilitation after surgery
- Protection during and after sports activities
- In rehabilitation after plaster cast

ONE
SIZE

COLOR ☒ BLACK ☐ WHITE



364 Padded Ankle Support

It is a supportive and comfortable ankle brace with a sponge pad designed in accordance with the joint.

Indications

- In recurrent ankle sprains in ligament damage
- Rehabilitation after surgery
- Protection during and after sports activities
- In rehabilitation after plaster cast

ONE
SIZE

COLOR ☒ BLACK ☐ WHITE



365 R.O.M WALKER

The base part is made of hard plastic material and the outer part is covered with a pile pad. It is an ankle brace that helps the patient's gait.

Indications

- It is used for postoperative stabilization in tendon, ligament and soft tissue injuries
- After acid tendon rupture surgery
- After heel and malleolar bone fractures



S	M	L
35-37	38-40	40-44
PERIMETER/ CM		

ANKLE



366 ACHILLES BOOT

It is an ankle brace orthosis that is made of EVA on the sole and hard plastic on the outside, covered with a terrycloth. It has a heel wedge and a strap across the foot to provide support and stability.

Indications

- Chronic ligament instability
- Tendon, ligament and soft tissue injuries
- Achilles tendon tears for postoperative stabilization
- After surgery or injury



S	M	L
34-37	38-41	42-45
PERIMETER/ CM		



367 Rest Afo (Foot Guard)

It is made of thermoplastic material and is designed in accordance with the foot anatomy. It is a T-band orthosis that provides support and stability to the foot.

Indications

- In chronic fasciitis syndrome
- Maintaining the foot position of the inpatient
- It is also used in rehabilitation after surgical intervention



S	M	L
35-37	38-40	40-44
PERIMETER/ CM		



368 Foot Guard Lux (Degree Adjustable)

It is made of thermoplastic material and is designed in accordance with the foot anatomy. It is a T-band orthosis with a heel wedge and a strap across the foot. The surface is covered with sponge fabric. Can be used for right and left feet. The non-slip surface on the sole provides ease of use for short distance walks.

Indications

- Plantar fasciitis
- Maintaining the position of the foot in the inpatient



ADJUSTABLE		
S	M	L
35-37	38-40	40-44
PERIMETER/ CM		



369 Reflex Afo (Right/Left)

It is an ankle orthosis made of plastic material and designed in accordance with the foot anatomy. It increases lateral stability and provides support and stability to the foot.

Indications

- Low foot deformity
- Muscle weakness in paralyzed patients
- Rehabilitation after surgery
- It is used in cases where it is desired to keep the ankle at 90 degrees



S	M	L
34-37	38-41	42-45
PERIMETER/ CM		



421 Hallux Valgus Splint

-t Ms a hallY\ ZaKYs sTIMRt HeZeloTeH to TreZeRt the TroKressMoR oJ hallY\ ZaKYs HeJorQMt]. -t coRsMsts oJ a thumbguardandelastic bandagescovering thethumb.Itprevents lateraldeviationbyapplyingapulling force to the thYQF.)lastMc aRH soJt QaterMal strYctYre Hoes Rot Jeel YRcoQJortaFle aRH TroZMHes Qa\M mumcomfort.Suitable forshoesanddaily use.

Indications

g ,allY\ ZaKYs HeJorQMt]

ONE
SIZE



422 Hallux Valgus Night Splint

Itismadeofspongefabricandvelcrobandsprovideeasyusageandstabilization.There is a support at the hallY\ ZaKYs area TroZMHes correct TosMtMoR Jor toe. -t Ms YseH MR ZaKYs HeJorQMt].

Indications

g ,allY\ ZaKYs HeJorQMt]

•Rehabilitationaftersurgical intervention

ONE
SIZE



423 Hallux Valgus Splint PLASTIC

-t Ms a hallY\ ZaKYs sYTTor t that aHNYsts 5e\MoR aRH e\tractMoR oJ the FMK toe. -t caR Fe easMI] attacheH aRH re-QoZeH aRH caR Fe YseH [Mth shoes. 8he thYQF Ms eJjectMZel] sYTTorteH. ,ealMRK ajter oTeratMoR Ms accelerateH.

Indications

AÇIKLAMA GELECEK

ONE
SIZE




FOOT HEALTH



400

SILICONeInsole



401

ArcHighSilicone Insole



402

SiliconeInsole Fabric Coated



403

MedialandLateral Wedge



405

SiliconeHeel Cup



406

Silicone Heel Cup Fabric Coated with Adhesive



407

SiliconeToe Separator



408

SiliconeArch Support with Adhesive



409

SiliconeBunion Protector




FOOT HEALTH



411

Silicone Metatarsal Support



414

Silicone Stress Ball



416

Starflex Insole



417

CORK INSOLE



418

Cork Insole with Epin



419

STEEL INSOLE



424

Hallux Valgus Shoe Sağ-sol



425

PLASTER TOE



426

SABO SHOE

Compression Stockings

POINT
NAVIGATING YOUR PATH TO OPTIMAL HEALTH



DESCRIPTIONS

Varicose Socks manufactured from quality yarns and they are produced according to world standards and common sizes. It provides maximum comfort due to technological knitting technique and yarn quality and it provides maximum comfort and also it enables to minimize sweating to prevent bacterial growth thanks to respirable feature. Drop silicones, prevent slipping and does not disturb. Materials used in knitting are harmless to body and environment



551-555

CLASS I - CLASS II

VARICOSE STOCKING ABOVE KNEE

COLOR ● BLACK ● BEIGE



552-554

CLASS I - CLASS II

VARICOSE STOCKING PANTYHOSE

COLOR ● BLACK ● BEIGE



556

CLASS I - CLASS II

VARICOSE STOCKING PREGNANCY

COLOR ● BLACK ● BEIGE



550-553

CLASS I - CLASS II

VARICOSE STOCKING BELOW KNEE

COLOR ● BLACK ● BEIGE

Varicose
STOCKINGS



Leg Circumference	SIZES						
	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cF	36-46	37-51	42-57	45-60	49-65	52-68	55-70
cG	43-58	46-61	49-64	52-68	46-71	60-74	64-79

Antiembolism STOCKINGS





557
ANTI-EMBOLISM STOCKING
BELOW KNEE



558
ANTI EMBOLISM STOCKING
ABOVE KNEE



559
ANTI EMBOLISM STOCKING
BELTED



Leg Circumference	SIZES						
	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cF	36-46	37-51	42-57	45-60	49-65	52-68	55-70
cG	43-58	46-61	49-64	52-68	46-71	60-74	64-79



450

SILICONE BREAST PROTHESIS



281

NoseGuard Mask



371

Suspensory Slip



373

Frejka Pillow



375

Pavlic Harness



452

TOILET SEAT RISER



700

Chamber Inhaler



374

DKC ORTHOSIS




376

DENNISBROWN SPLINT



MOTIIA





600

VISCO NECK PILLOW



601

VISCONECKGEL PILLOW



602

VISCOTRAVEL PILLOW



603

ViscoRingSeat Cushion
607-FOAM Ring Seat Cushion




604

Visco Back Cushion
609-FOAM visco back cushion



605

ViscoAnatomic Seat Cushion



606

REFLUX PILLOW



608

RingSeatCushion (OPEN)



610

Position Pillow



NAVIGATING YOUR PATH TO OPTIMAL HEALTH

WALKING AID



WALKING AID




650

WALKING STICK



651

TRIPODWALKING STICK




653

CRUTCHES



654

CRUTCHES LUX



655

CANEDIAN CRUTCHES



656

WALKER (IRON)



657

WALKER (ALUMINUM)



658

FOLDINGCOMODE



661

MANUAL WHEELCHAIR



